BAUSCH+ Health

DESCRIPTION: E5699 E5699P Corneal Protector Ins / GBL / STL

PART No.: 4119903

SPEC No. or SPEC DIMENSIONS: REF STL p/n S21507

SPECIAL INSTRUCTIONS: n/a

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CROUCH CORNEAL PROTECTOR E5699 & E5699 P

ENG

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Description:

The Crouch Corneal Protector has been developed for oculo-plastic work. The overall dimensions are as follows:

	Adult	Pediatric
Horizontal inner diameter	24.2mm	16.0mm
Horizontal outer diameter	25.8mm	17.0mm
Vertical inside diameter	22.0mm	16.0mm
Vertical outside diameter	23.5mm	17.0mm
Inside Height	8.7mm	7.0mm

These dimensions enable the protector to rest on the sclera and not on the cornea. During plastic surgery, the cornea and sclera are protected from manipulation.

Indications:

Corneal and scleral protection is essential when performing oculo-plastic procedures. The Crouch Corneal Protector is indicated in the following conditions for surgical procedures:

1. Adult Blepharoplasty
2. Levator resections
3. Fascia lata brow suspension procedures
4. Orbital trauma
5. Entropion repair or ectropion repair of the lids and other plastic surgical procedures for the lid.

The corneal protector should be removed after surgery.

Contraindications:

1. Severe microphthalmus
2. Definite or suspected rupture of the globe.
3. Corneal laceration.

MEDICAL DEVICE RE-USE STATEMENT

If this product is reprocessed and/or re-used, Bausch + Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient or user and, in extreme incidents, death. This product is labeled as 'single use' which is defined as a device intended to be used once only for a single patient.

Instructions for Use:

Normally a lubricant is not required for this device. However, if desired, the anterior and posterior surfaces of the protector may be covered with a thin layer of a lubricating agent.

The corneal protector is inserted by first grasping the tab on the protector and sliding it under the superior lid. The lower lid is then retracted and the corneal protector is placed into position on the sclera. The corneal protector can be easily removed after surgery with a Stevens Hook (E0600) placed beneath the protector. If a muscle hook is not available, the corneal protector can be removed by grasping the tab on the protector with a Bishop Harmon Forceps (E1500).

It is essential to remove the corneal protector after surgery.

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